

October 7, 2009
DRAFT

INTEGRATING THE HEALTHCARE ENTERPRISE (IHE) PRE-CONNECTATHON TESTING

Patient Identity Cross Reference (PIX) and Patient Discovery Query (PDQ) Test Tool User's Guide Supplement

Submitted by:

Rob Snelick
National Institute of Standards and Technology (NIST)

TABLE OF CONTENTS

1.0	OVERVIEW	3
2.0	SUPPORT FILES	4
2.1	HL7 V2 Messages	4
2.1.1	Test Case 10501: PIX Query Case 1 and Case 2	4
2.1.2	Test Case 10502: PIX Query Case 3.....	11
2.1.3	Test Case 10503: PIX Query Case 4.....	18
2.1.4	Test Case 10506: PIX Query, Patient Update	23
2.1.5	Test Case 10511: PIX Patient Feed: A01	32
2.1.6	Test Case 10512: PIX Patient Feed: A04	40
2.1.7	Test Case 10515: PIX Patient Feed: A40	42
2.2	HL7 V3 Messages	48
2.2.1	Test Case 10501v3: PIX Query Case 1 and Case 2	48
2.2.2	Test Case 10502v3: PIX Query Case 3	58
2.2.3	Test Case 10503v3: PIX Query Case 4.....	65
2.2.4	Test Case 10506v3: PIX Query Patient Update	67
2.2.5	Test Case 10512v3: PIX Patient Feed.....	76
2.2.6	Test Case 10515v3: Patient Identity Merge	81

1.0 OVERVIEW

This is a supplement to the NIST PIX/PDQ Pre-Connectation Test Tool. This document contains test messages and the validation context files. The validation context files contain the list of assestions not captured the general conformance testing mechanisms.

2.0 SUPPORT FILES

2.1 HL7 V2 Messages

2.1.1 TEST CASE 10501: PIX QUERY CASE 1 AND CASE 2

Message Type	Message
Acknowledge Patient Add: Number 1	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20090810140837 ACK^A04^ACK NIST-090810140837273 P 2.3.1 MSA AA Flame-990001
Acknowledge Patient Add: Number 2	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20090810140842 ACK^A04^ACK NIST-090810140842899 P 2.3.1 MSA AA NIST-090810140838539
Acknowledge Patient Add: Number 3	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20090810140848 ACK^A04^ACK NIST-090810140848555 P 2.3.1 MSA AA NIST-090810140844164
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> <n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> <FailureInterpretation> <MessageFailure Type="VERSION" Level="ERROR" /> <MessageFailure Type="LENGTH" Level="IGNORE" /> <MessageFailure Type="DATATYPE" Level="ERROR" /> <MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> <MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> <MessageFailure Type="DATA" Level="ERROR" /> <MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> <MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

Message Type	Message
	<pre> - <!-- MSA.1 --> = <DataValueLocationItem> = <Location> = <Segment Name="MSA"> <Field Position="1" /> </Segment> </Location> = <Value> <PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 2009081 0141748 ADT^A04^ADT_A01 Flame-990001 P 2.3.1 EVN A04 20090810 PID PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI ALPHA^AL AN Barnes 19781208 M 820 JORIE BLVD^^New York City^NY^10503 (112)555-1234 (112)555-6789 153-12-5432 PV1 V </pre>
Patient Add: Number 2	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 2009081 0141755 ADT^A04^ADT_A01 NIST-090810141751293 1 2.3.1 EVN A04 20090810 PID XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI ALPHA^ALA N Barnes 19781208 M 820 JORIE BLVD^^New York City^NY^10503 (112)555-1234 (112)555-6789 153-12-5432 PV1 V </pre>
Patient Add: Number 3	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 2009081 0141800 ADT^A04^ADT_A01 NIST-090810141756903 1 2.3.1 EVN A04 20090810 PID ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI SIMPSON^ CARL AMSTRONG 19781209 M 820 OREL BLVD^^CHICAGO^IL^60523 (312)555- 1234 (312)555-6789 163-12-5432 PV1 V </pre>
PIX Query: Number 1	<pre> MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST 20090 730114101 QBP^Q23^QBP_Q21 NIST-090730114100989 D 2.5 QPD IHE PIX QUERY QRY1248968460880 PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1 &ISO^PI XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO RCP I </pre>
PIX Query: Number 2	<pre> MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST 20090 526133934 QBP^Q23^QBP_Q21 NIST-090526133933407 D 2.5 QPD IHE PIX QUERY QRY1243359573282 ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1 &ISO^PI text^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO RCP I </pre>
PIX Query Response: Number	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST 20090 728151928 RSP^K23^RSP_K23 NIST- </pre>

Message Type	Message
1	090728151927498 P 2.5 124.3.11.62382.03 MSA AA NIST-090728151926716 QAK QRY1248808766638 OK QPD IHE PIX QUERY QRY1248808766638 PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1 &ISO XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO PID XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI ~^^^S
PIX Query Response: Number 2	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST 20090 526133935 RSP^K23^RSP_K23 NIST- 090526133935407 P 2.5 124.3.11.62382.03 MSA AA NIST-090526133933407 QAK QRY1243359573282 NF QPD IHE PIX QUERY QRY1243359573282 ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1 &ISO^PI text^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO
PIX Query Response Validation Context: Number 1	<?xml version="1.0" encoding="UTF-8" ?> = <n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/message" > = <FailureInterpretation> <MessageFailure Type="VERSION" Level="ERROR" /> <MessageFailure Type="LENGTH" Level="IGNORE" /> <MessageFailure Type="DATATYPE" Level="ERROR" /> <MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> <MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> <MessageFailure Type="DATA" Level="ERROR" /> <MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> = <MessageInstanceSpecificValues> - <!-- PIX Manager RSP K23 --> - <!-- MSH 9.2 --> = <DataValueLocationItem> = <Location> = <Segment Name="MSH"> = <Field Position="9"> <Component Position="2" /> </Field> </Segment> </Location> = <Value> <PlainText>K23</PlainText> </Value> </DataValueLocationItem> - <!-- MSH 9.3

Message Type	Message
	<pre> --> = <DataValueLocationItem> = <Location> = <Segment Name="MSH"> = <Field Position="9"> = <Component Position="3" /> = </Field> = </Segment> = </Location> = <Value> = <PlainText>RSP_K23</PlainText> = </Value> = </DataValueLocationItem> - <!-- PID 3.1 --> = <DataValueLocationItem> = <Location> = <Segment Name="PID"> = <Field Position="3"> = <Component Position="1" /> = </Field> = </Segment> = </Location> = <Value> = <PlainText>XYZ10501</PlainText> = </Value> = </DataValueLocationItem> - <!-- PID 3.4.1 --> = <DataValueLocationItem> = <Location> = <Segment Name="PID"> = <Field Position="3"> = <Component Position="4"> = <SubComponent Position="1" /> = </Component> = </Field> = </Segment> = </Location> = <Value> = <PlainText>XREF2005</PlainText> = </Value> = </DataValueLocationItem> - <!-- PID 3.4.2 --> = <DataValueLocationItem> </pre>

Message Type	Message
	<pre> = <Location> = <Segment Name="PID"> = <Field Position="3"> = <Component Position="4"> = <SubComponent Position="2" /> </Component> </Field> </Segment> </Location> = <Value> = <PlainText>1.3.6.1.4.1.21367.2005.1.2</PlainText> </Value> </DataValueLocationItem> - <!-- PID 3.4.3 --> = <DataValueLocationItem> = <Location> = <Segment Name="PID"> = <Field Position="3"> = <Component Position="4"> = <SubComponent Position="3" /> </Component> </Field> </Segment> </Location> = <Value> = <PlainText>ISO</PlainText> </Value> </DataValueLocationItem> - <!-- PID 5.7 --> = <DataValueLocationItem> = <Location> = <Segment Name="PID"> = <Field Position="5"> = <Component Position="7" /> </Field> </Segment> </Location> = <Value> = <PlainText>S</PlainText> </Value> </DataValueLocationItem> - <!-- MSA 1 --> = <DataValueLocationItem> </pre>

Message Type	Message
	<pre> => <Location> => <Segment Name="MSA"> => <Field Position="1" /> => </Segment> => </Location> => <Value> => <PlainText>AA</PlainText> => </Value> => </DataValueLocationItem> - <!-- QAK 2 --> => <DataValueLocationItem> => <Location> => <Segment Name="QAK"> => <Field Position="2" /> => </Segment> => </Location> => <Value> => <PlainText>OK</PlainText> => </Value> => </DataValueLocationItem> => </MessageInstanceSpecificValues> => </n1:HL7V2MessageValidationContextDefinition> </pre>
PIX Query Response Validation Context: Number 2	<pre> <?xml version="1.0" encoding="UTF-8" ?> => <n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/message" > => <FailureInterpretation> => <MessageFailure Type="VERSION" Level="ERROR" /> => <MessageFailure Type="LENGTH" Level="IGNORE" /> => <MessageFailure Type="DATATYPE" Level="ERROR" /> => <MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> => <MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> => <MessageFailure Type="DATA" Level="ERROR" /> => <MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> => </FailureInterpretation> => <MessageInstanceSpecificValues> - <!-- PIX Manager RSP K23 --> - <!-- MSH 9.1 --> => <DataValueLocationItem> => <Location> => <Segment Name="MSH"> => <Field Position="9"> </pre>

Message Type	Message
	<pre> <Component Position="1" /> </Field> </Segment> </Location> = <Value> <PlainText>RSP</PlainText> </Value> </DataValueLocationItem> - <!-- MSH 9.2 --> = <DataValueLocationItem> = <Location> = <Segment Name="MSH"> = <Field Position="9"> = <Component Position="2" /> </Field> </Segment> </Location> = <Value> <PlainText>K23</PlainText> </Value> </DataValueLocationItem> - <!-- MSH 9.3 --> = <DataValueLocationItem> = <Location> = <Segment Name="MSH"> = <Field Position="9"> <Component Position="3" /> </Field> </Segment> </Location> = <Value> <PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> - <!-- MSA 1 --> = <DataValueLocationItem> = <Location> = <Segment Name="MSA"> <Field Position="1" /> </Segment> </Location> = <Value> <PlainText>AA</PlainText> </pre>

Message Type	Message
	<pre> </Value> </DataValueLocationItem> - <!-- QAK 2 --> = <DataValueLocationItem> = <Location> = <Segment Name="QAK"> = <Field Position="2" /> </Segment> </Location> = <Value> <PlainText>NF</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

2.1.2 TEST CASE 10502: PIX QUERY CASE 3

Message Type	Message
Acknowledge Patient Add: Number 1	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942731650 ACK^A04^ACK 12434334316500 P 2.3.1 MSA AA NIST-090527101025712 </pre>
Acknowledge Patient Add: Number 2	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942736416 ACK^A04^ACK 12434334364161 P 2.3.1 MSA AA NIST-090527101031978 </pre>
Acknowledge Validation Content: All	<pre> : ?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- </pre>

Message Type	Message
	<pre> <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527101029 ADT^A04^ADT_A01 NIST-090527101025712 1 2.3.1 EVN A04 20090527 PID PIX10502^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI B ETA^BETTY BRAYAN 19781208 F 820 JORIE BLVD^^CHICAGO^IL^60523 (312)555-1234 (312)555-6789 163-12- 5432 PV1 V </pre>
Patient Add: Number 2	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527101035 ADT^A04^ADT_A01 NIST-090527101031978 1 2.3.1 EVN A04 20090527 PID XYZ10502^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI C ROSS^KEN SMITH 19781209 M 820 OREL BLVD^^CHICAGO^IL^60523 (312)555-1234 (312)555-6789 163-12- 5432 PV1 V </pre>
PIX Query	<pre> MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIS T 20090527102956 QBP^Q23^QBP_Q21 NIST-090527102956719 D 2.5 QPD IHE PIX QUERY QRY1243434596641 XX10502^^^HIMSS2005&1.3.6.1.4.1.21367.20 05.1.1&ISO RCP I </pre>
PIX Query	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIS </pre>

Message Type	Message
Response	T 20090527102958 RSP^K23^RSP_K23 NIST- 090527102958219 P 2.5 124.3.11.62382.03 MSA AE NIST-090527102956719 Invalid patient ERR QPD^1^3 QPD^1^3^1^1 204 W QAK QRY1243434596641 AE QPD IHE PIX QUERY QRY1243434596641 XX10502^^^HIMSS2005&1.3.6.1.4.1.21367.20 05.1.1&ISO
PIX Query Response Validation Context	<pre> : ?xml version="1.0" encoding="UTF-8" ?> n1: HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> </pre>

Message Type	Message
	<pre> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>AE</PlainText> </pre>

Message Type	Message
	<pre> </Value> </DataValueLocationItem> <!-- ERR.1.1 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="1"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>QPD</PlainText> </Value> </DataValueLocationItem> <!-- ERR.1.2 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="1"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR.1.3 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="1"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>3</PlainText> </Value> </DataValueLocationItem> <!-- </pre>

Message Type	Message
	<pre> ERR 2.1 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>QPD</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.2 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.3 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>3</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.4 --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> Location> Segment Name="ERR"> Field Position="2"> Component Position="4" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.5 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="5" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR 3.1 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>204</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

2.1.3 TEST CASE 10503: PIX QUERY CASE 4

Message Type	Message
PIX Query Case	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST T 20090527111441 QBP^Q23^QBP_Q21 NIST-090527111440914 D 2.5 QPD IHE PIX QUERY QRY1243437280836 ABC10503^^^XXXX&1.3.6.1.4.1.21367.2005.3 .3333&ISO RCP I
PIX Query Response	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST T 20090527111442 RSP^K23^RSP_K23 NIST- 090527111442320 P 2.5 124.3.11.62382.03 MSA AE NIST-090527111440914 Invalid assigning authority ERR QPD^1^3 QPD^1^3^1^4 204 W QAK QRY1243437280836 AE QPD IHE PIX QUERY QRY1243437280836 ABC10503^^^XXXX&1.3.6.1.4.1.21367.2005.3 .3333&ISO
PIX Query Response Validation Context	<?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value>

Message Type	Message
	<pre> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- QPD 3.1 --> DataValueLocationItem> Location> Segment Name="QPD"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>ABC10503</PlainText> </Value> </DataValueLocationItem> <!-- QPD 3.4.1 </pre>

Message Type	Message
	<pre> --> DataValueLocationItem> Location> Segment Name="QPD"> Field Position="3"> Component Position="4"> SubComponent Position="1" /> </Component> </Field> </Segment> </Location> Value> PlainText>XXXX</PlainText> </Value> </DataValueLocationItem> <!-- QPD 3.4.2 --> DataValueLocationItem> Location> Segment Name="QPD"> Field Position="3"> Component Position="4"> SubComponent Position="2" /> </Component> </Field> </Segment> </Location> Value> PlainText>1.3.6.1.4.1.21367.2005.3.3333</PlainText> </Value> </DataValueLocationItem> <!-- QPD 3.4.3 --> DataValueLocationItem> Location> Segment Name="QPD"> Field Position="3"> Component Position="4"> SubComponent Position="3" /> </Component> </Field> </Segment> </Location> Value> PlainText>ISO</PlainText> </pre>

Message Type	Message
	<pre> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.1 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>QPD</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.2 --> DataValueLocationItem> Location> </pre>

Message Type	Message
	<pre> Segment Name="ERR"> Field Position="2"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.3 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>3</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.4 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> Component Position="4" /> </Field> </Segment> </Location> Value> PlainText>1</PlainText> </Value> </DataValueLocationItem> <!-- ERR 2.5 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="2"> </pre>

Message Type	Message
	<pre> Component Position="5" /> </Field> </Segment> </Location> Value> PlainText>4</PlainText> </Value> </DataValueLocationItem> <!-- ERR 3.1 --> DataValueLocationItem> Location> Segment Name="ERR"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>204</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

2.1.4 TEST CASE 10506: PIX QUERY, PATIENT UPDATE

Message Type	Message
Acknowledge Patient Add: Number 1	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942737167 ACK^A04^ACK 12434381771670 P 2.3.1 MSA AA NIST-090527112931073 </pre>
Acknowledge Patient Add: Number 2	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20094274258 ACK^A04^ACK 12434381820581 P 2.3.1 MSA AA NIST-090527112937401 </pre>
Acknowledge Patient Update	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942723273 ACK^A08^ACK 12434466232732 P 2.3.1 MSA AA NIST-090527135016929 </pre>
Acknowledge Patient Update: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> </pre>

Message Type	Message
	<pre> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> </pre>

Message Type	Message
	<pre> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527112934 ADT^A04^ADT_A01 NIST-090527112931073 1 2.3.1 EVN A04 20090527 PID PIX10506^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI T AU^TERI BEVERE 06041555 M 202 KEN HABOR^^NEW YORK CITY^NY^61000 (213)555-1234 (213)500-1234 361-21-2345 PV1 V </pre>
Patient Add: Number 2	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527112941 ADT^A04^ADT_A01 NIST-090527112937401 1 2.3.1 EVN A04 20090527 PID ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI TO W^T BEVERE 19781555 M 022 HELEN GAD^^CHICAGO^IL^11000 (213)555-1234 (213)005-1234 113-005- 5678 PV1 V </pre>

Message Type	Message
Patient Update	MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527135020 ADT^A08^ADT_A01 NIST-090527135016929 1 2.3.1 EVN A08 20090527 PID ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI TA U^TERI BEVERE 06041555 M 202 KEN HAVOR^^NEW YORK CITY^NY^61000 (213)555-1234 (213)500-1234 361-21-2345 PV1 I
PIX Query: Number 1	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIS T 20090527113947 QBP^Q23^QBP_Q21 NIST-090527113946959 D 2.5 QPD IHE PIX QUERY QRY1243438786881 PIX10506^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&I SO RCP I
PIX Query: Number 2	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIS T 20090527135721 QBP^Q23^QBP_Q21 NIST-090527135721677 D 2.5 QPD IHE PIX QUERY QRY1243447041583 PIX10506^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&I SO RCP I
PIX Query Response: Number 1	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIS T 20090527113948 RSP^K23^RSP_K23 NIST- 090527113948444 P 2.5 124.3.11.62382.03 MSA AA NIST-090527113946959 QAK QRY1243438786881 NF QPD IHE PIX QUERY QRY1243438786881 PIX10506^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&I SO
PIX Query Response: Number 2	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIS T 20090527135723 RSP^K23^RSP_K23 NIST- 090527135723286 P 2.5 124.3.11.62382.03 MSA AA NIST-090527135721677 QAK QRY1243447041583 OK QPD IHE PIX QUERY QRY1243447041583 PIX10506^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&I SO PID ABC10506^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI ~^ ^^^^S
PIX Query Response Validation Context: Number 1	<?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" />

Message Type	Message
	<pre> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </pre>

Message Type	Message
	<pre> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>NF</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
PIX Query Response Validation Context: Number 2	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> </pre>

Message Type	Message
	<pre> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </pre>

Message Type	Message
	<pre> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>OK</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.1 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>ABC10506</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.1 --> DataValueLocationItem> Location> Segment Name="PID"> </pre>

Message Type	Message
	<pre> Field Position="3"> Component Position="4"> SubComponent Position="1" /> </Component> </Field> </Segment> </Location> Value> PlainText>XREF2005</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.2 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="2" /> </Component> </Field> </Segment> </Location> Value> PlainText>1.3.6.1.4.1.21367.2005.1.2</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.3 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="3" /> </Component> </Field> </Segment> </Location> Value> PlainText>ISO</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.5 </pre>

Message Type	Message
	<pre>--> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="5" /> </Field> </Segment> </Location> Value> PlainText>PI</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition></pre>

2.1.5 TEST CASE 10511: PIX PATIENT FEED: A01

Message Type	Message
Acknowledge Patient Add: Number 1	<pre>MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20094275105 ACK^A01^ACK 12434482251050 P 2.3.1 MSA AA Flame-990001 </pre>
Acknowledge Patient Add: Number 2	<pre>MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 20094279652 ACK^A01^ACK 12434482296521 P 2.3.1 MSA AA NIST-090527141705324 </pre>
Acknowledge Patient Add: Number 3	<pre>MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942714996 ACK^A01^ACK 12434482349962 P 2.3.1 MSA AA NIST-090527141710715 </pre>
Acknowledge Validation Content: All	<pre><?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --></pre>

Message Type	Message
	<pre> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527141702 ADT^A01^ADT_A01 Flame-990001 P 2.3.1 EVN A01 20090527 PID PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI A LPHA^ALAN Barnes 19781208 M 820 JORIE BLVD^^New York City^NY^10503 (112)555-1234 (112)555-6789 153-12-5432 PV1 V </pre>
Patient Add: Number 2	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527141708 ADT^A01^ADT_A01 NIST-090527141705324 A 2.3.1 EVN A01 20090527 PID XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI AL PHA^ALAN Barnes 19781208 M 820 JORIE BLVD^^New York City^NY^10503 (112)555-1234 (112)555-6789 153-12-5432 PV1 V </pre>
Patient Add: Number 3	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527141714 ADT^A01^ADT_A01 NIST-090527141710715 A 2.3.1 EVN A01 20090527 PID ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI S IMPSON^CARL AMSTRONG 19781209 M 820 OREL BLVD^^CHICAGO^IL^60523 (312)555-1234 (312)555-6789 163-12- 5432 PV1 V </pre>

Message Type	Message
PIX Query: Number 1	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST T 20090527142638 QBP^Q23^QBP_Q21 NIST-090527142637983 D 2.5 QPD IHE PIX QUERY QRY1243448797905 PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO RCP I
PIX Query: Number 2	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST T 20090527142639 QBP^Q23^QBP_Q21 NIST-090527142639342 D 2.5 QPD IHE PIX QUERY QRY1243448799280 ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI text^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO RCP I
PIX Query Response: Number 1	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST T 20090527142639 RSP^K23^RSP_K23 NIST- 090527142638655 P 2.5 124.3.11.62382.03 MSA AA NIST-090527142637983 QAK QRY1243448797905 OK QPD IHE PIX QUERY QRY1243448797905 PIX10501^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO PID XYZ10501^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI ~^ ^^^^^S
PIX Query Response: Number 2	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST T 20090527142641 RSP^K23^RSP_K23 NIST- 090527142641077 P 2.5 124.3.11.62382.03 MSA AA NIST-090527142639342 QAK QRY1243448799280 NF QPD IHE PIX QUERY QRY1243448799280 ABC10501^^^HIMSS2005&1.3.6.1.4.1.21367.2 005.1.1&ISO^PI text^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO
PIX Query Response Validation Context: Number 1	<?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1

Message Type	Message
	<pre> --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.1 --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> Location> Segment Name="PID"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>XYZ10501</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.1 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="1" /> </Component> </Field> </Segment> </Location> Value> PlainText>XREF2005</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.2 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="2" /> </Component> </Field> </Segment> </Location> Value> PlainText>1.3.6.1.4.1.21367.2005.1.2</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.3 </pre>

Message Type	Message
	<pre> --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="3" /> </Component> </Field> </Segment> </Location> Value> PlainText>ISO</PlainText> </Value> </DataValueLocationItem> <!-- PID 5.7 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="5"> Component Position="7" /> </Field> </Segment> </Location> Value> PlainText>S</PlainText> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>OK</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
PIX Query Response Validation Context: Number 2	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.2 --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </pre>

Message Type	Message
	<pre> </Location> Value> PlainText>NF</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

2.1.6 TEST CASE 10512: PIX PATIENT FEED: A04

Message Type	Message
Acknowledge Patient Add: Number 1	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942735276 A CK^A04^ACK 12434560552760 P 2.3.1 MSA AA NIST-090527162729385 </pre>
Acknowledge Patient Add: Number 2	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942735276 A CK^A04^ACK 12434560552760 P 2.3.1 MSA AA NIST-090527162729385 </pre>
Acknowledge Patient Add: Number 3	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942735276 A CK^A04^ACK 12434560552760 P 2.3.1 MSA AA NIST-090527162729385 </pre>
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </pre>

Message Type	Message
	<pre> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527144212 ADT^A04^ADT_A01 NIST-090527144209435 1 2.3.1 EVN A04 20090527 PID PIX10512^^^HIMSS2005 EPSILON^ELLIE JACK 19380224 F 1 PINETREE^^WEBSTER^MO^63119 (314)555-1234 (314)555- 4444 163-12-5432 PV1 V </pre>
Patient Add: Number 2	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527144846 ADT^A04^ADT_A01 NIST-090527144842416 1 2.3.1 EVN A04 20090527 PID PIX10512^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI E PSILON^ELLIE JACK 19380224 F 1 PINETREE^^WEBSTER^MO^63119 (314)555-1234 (314)555- 4444 163-12-5432 PV1 V </pre>
Patient Add: Number 3	<pre> MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090527162732 ADT^A04^ADT_A01 NIST-090527162729385 1 2.3.1 EVN A04 20090527 PID PIX10512^^^&1.3.6.1.4.1.21367.2005.1.1&ISO EPSILON^ELLIE JACK 19380224 F 1 PINETREE^^WEBSTER^MO^63119 (314)555- 1234 (314)555-4444 163-12-5432PV1 V </pre>

2.1.7 TEST CASE 10515: PIX PATIENT FEED: A40

Message Type	Message
Acknowledge Patient Add: Number 1	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942828226 ACK^A04^ACK 12435227682260 P 2.3.1 MSA AA NIST-090528105921882 </pre>
Acknowledge	<pre> MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST </pre>

Message Type	Message
Patient Add: Number 2	200942833507 ACK^A04^ACK 12435227735071 P 2.3.1 MSA AA NIST-090528105928523
Acknowledge Patient Add: Number 3	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942839289 ACK^A04^ACK 12435227792892 P 2.3.1 MSA AA NIST-090528105934554
Acknowledge Patient Merge	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Flame_PIX_Source NIST 200942826665 ACK^A40^ACK 12435228266653 P 2.3.1 MSA AA NIST-090528110022806
Acknowledge Content Validation: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- MSH.11.1 --> <!-- <DataValueLocationItem> <Location> <Segment Name="MSH"> <Field Position="11"> <Component Position="1"> </Component> </Field> </Segment> </Location> <Value><PlainText>P</PlainText></Value> </DataValueLocationItem> --> <!-- MSA.1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> </pre>

Message Type	Message
	PlainText>AA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition>
Patient Add: Number 1	MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090528105925 ADT^A04^ADT_A01 NIST-090528105921882 P 2.3.1 EVN A04 20090528 PID PIX10515W^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI WASHINGTON^MARY JOHN 19781208 F 820 JORIE BLVD^^CHICAGO^IL^60523 (312)555-1234 (312)555-6789 100-9-1234 PV1 v
Patient Add: Number 2	MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090528105932 ADT^A04^ADT_A01 NIST-090528105928523 P 2.3.1 EVN A04 20090528 PID XYZ10515W^^^XREF2005&1.3.6.1.4.1.21367.2005.1.2&ISO^PI W ASHINGTON^MARY JOHN 19781208 F 820 JORIE BLVD^^CHICAGO^IL^60523 (312)555-1234 (312)555-6789 100-9-1234 PV1 v
Patient Add: Number 3	MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090528105938 ADT^A04^ADT_A01 NIST-090528105934554 P 2.3.1 EVN A04 20090528 PID PIX10515L^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI LINCOLN^MARY JOHN 19781208 F JEAN JAQUES BLVD^^New York^NY^60001 (213)455-4321 (213)555-0001 100-9-1234 PV1 v
Patient Merge: Number 1	MSH ^~\& NIST_Flame_PIX_Source NIST NIST_Swan_PIX_Manager NIST 20090528110025 ADT^A40^ADT_A39 NIST-090528110022806 P 2.3.1 EVN A40 20090528 PID PIX10515L^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI LINCOLN^MARY JOHN 19781208 F JEAN JAQUES BLVD^^New York^NY^60001 (213)455-4321 (213)555-0001 100-9-1234 MRG PIX10515W^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO PID ^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI LINCOLN^M ARY JOHN 19781208 F JEAN JAQUES BLVD^^New York^NY^60001 (213)455-4321 (213)555-0001 MRG PIX10515W^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO PID ^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI LINCOLN^M ARY JOHN 19781208 F JEAN JAQUES BLVD^^New York^NY^60001 (213)455-4321 (213)555-0001 MRG PIX10515W^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO PID ^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI LINCOLN^M ARY JOHN 19781208 F JEAN JAQUES BLVD^^New York^NY^60001 (213)455-4321 (213)555-0001 MRG PIX10515W^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO

Message Type	Message
PIX Query	MSH ^~\& NIST_Hydra_PIX_Consumer NIST NIST_Swan_PIX_Manager NIST T 20090528110358 QBP^Q23^QBP_Q21 NIST-090528110358016 P 2.5 QPD IHE PIX QUERY QRY1243523037937 XYZ10515W^^^XREF2005&1.3.6.1.4.1.21367.2 005.1.2&ISO RCP I
PIX Query Response	MSH ^~\& NIST_Swan_PIX_Manager NIST NIST_Hydra_PIX_Consumer NIST T 20090528110359 RSP^K23^RSP_K23 NIST- 090528110359359 P 2.5 124.3.11.62382.03 MSA AA NIST-090528110358016 QAK QRY1243523037937 OK QPD IHE PIX QUERY QRY1243523037937 XYZ10515W^^^XREF2005&1.3.6.1.4.1.21367.2 005.1.2&ISO PID PIX10515L^^^HIMSS2005&1.3.6.1.4.1.21367.2005.1.1&ISO^PI ~^~~~~S
PIX Query Response Validation Context	<?xml version="1.0" encoding="UTF-8" ?> n1:HL7V2MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- PIX Manager RSP K23 --> <!-- MSH 9.1 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>RSP</PlainText> </Value> </DataValueLocationItem> <!--

Message Type	Message
	<pre> MSH 9.2 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="2" /> </Field> </Segment> </Location> Value> PlainText>K23</PlainText> </Value> </DataValueLocationItem> <!-- MSH 9.3 --> DataValueLocationItem> Location> Segment Name="MSH"> Field Position="9"> Component Position="3" /> </Field> </Segment> </Location> Value> PlainText>RSP_K23</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.1 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="1" /> </Field> </Segment> </Location> Value> PlainText>PIX10515L</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.1 --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="1" /> </Component> </Field> </Segment> </Location> Value> PlainText>HIMSS2005</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.2 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="2" /> </Component> </Field> </Segment> </Location> Value> PlainText>1.3.6.1.4.1.21367.2005.1.1</PlainText> </Value> </DataValueLocationItem> <!-- PID 3.4.3 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="3"> Component Position="4"> SubComponent Position="3" /> </Component> </Field> </Segment> </Location> Value> PlainText>ISO</PlainText> </Value> </DataValueLocationItem> </pre>

Message Type	Message
	<pre> <!-- PID 5.7 --> DataValueLocationItem> Location> Segment Name="PID"> Field Position="5"> Component Position="7" /> </Field> </Segment> </Location> Value> PlainText>S</PlainText> </Value> </DataValueLocationItem> <!-- MSA 1 --> DataValueLocationItem> Location> Segment Name="MSA"> Field Position="1" /> </Segment> </Location> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- QAK 2 --> DataValueLocationItem> Location> Segment Name="QAK"> Field Position="2" /> </Segment> </Location> Value> PlainText>OK</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V2MessageValidationContextDefinition> </pre>

2.2 HL7 V3 Messages

2.2.1 TEST CASE 10501V3: PIX QUERY CASE 1 AND CASE 2

Message Type	Message
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> DataValueLocationItem> <XPath>/*:MCCI_IN000002UV01/*:acknowledgement/*:typeCode /@code</XPath> Value> PlainText>CA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="b0b596f4-e145-4ad2-8e12-d20c623960af" /> creationTime value="20090818163443" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </pre>

Message Type	Message
	<pre> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="100" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>ALPHA</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Barnes</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> </pre>

Message Type	Message
	<pre> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Patient Add: Number 2	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="bf433029-255c-481c-8966-4e718cb3208c" /> creationTime value="20090818162938" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> </pre>

Message Type	Message
	<pre> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="200" assigningAuthorityName="domain2" root="1.2.3.4.5.2000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>ALPHA</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Barnes</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </pre>

Message Type	Message
	</PRPA_IN201301UV02>
Patient Add: Number 3	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:h17-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="02778280-e801-4ee2-94a8-d5eb79581d8c" /> creationTime value="20090818162949" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="101" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>SIMPSON</family> </name> </pre>

Message Type	Message
	<pre> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>AMSTRONG</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
PIX Query: Number 1	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="4a8aa794-b703-401c-8677-cca493c22d14" /> creationTime value="20090821104213" /> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" </pre>

Message Type	Message
	<pre> classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10501Q1" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> dataSource> value root="1.2.3.4.5.2000" /> semanticsText>DataSource.id</semanticsText> </dataSource> patientIdentifier> value extension="100" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </PRPA_IN201309UV02> </pre>
PIX Query: Number 2	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="b425af6a-0ade-4720-89fc-599c3bda76cf" /> creationTime value="20090821104228" /> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> </pre>

Message Type	Message
	<pre> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10501Q2" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> dataSource> value root="1.2.3.4.5.2000" /> semanticsText>DataSource.id</semanticsText> </dataSource> patientIdentifier> value extension="101" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </PRPA_IN201309UV02> </pre>
PIX Query Response Validation Context: Number 1	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> </pre>

Message Type	Message
	<pre> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- assigningAuthorityName --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@assigningAuthorit yName</XPath> Value> PlainText>domain2</PlainText> </Value> </DataValueLocationItem> <!-- extension --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@extension</XPath> > Value> PlainText>200</PlainText> </Value> </DataValueLocationItem> <!-- root --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@root</XPath> Value> PlainText>1.2.3.4.5.2000</PlainText> </Value> </DataValueLocationItem> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:typeCod e/@code</XPath> Value> PlainText>AA</PlainText> </pre>

Message Type	Message
	<pre> </Value> </DataValueLocationItem> <!-- Query ack code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryId/@root</XPath> Value> PlainText>10501Q1</PlainText> </Value> </DataValueLocationItem> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryResponseCode/@code</XPath> Value> PlainText>OK</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
PIX Query Response Validation Context: Number 2	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:typeCod e/@code</XPath> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- </pre>

Message Type	Message
	<pre> Query ack code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryId/@root</XPath> Value> PlainText>10501Q2</PlainText> </Value> </DataValueLocationItem> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryResponseCode/@code</XPath> Value> PlainText>NF</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>

2.2.2 TEST CASE 10502V3: PIX QUERY CASE 3

Message Type	Message
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> DataValueLocationItem> <XPath>/*:MCCI_IN000002UV01/*:acknowledgement/*:typeCode /@code</XPath> Value> </pre>

Message Type	Message
	<pre> PlainText>CA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="c3b31337-dd62-4b69-b70b-17bf67bf9a00" /> creationTime value="20090818163705" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="100" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> </patient> </subject1> </registrationEvent> </subject> </controlActProcess> </receiver> </interactionId> </PRPA_IN201301UV02> </pre>

Message Type	Message
	<pre> name> given>John</given> family>ALPHA</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Barnes</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Patient Add: Number 2	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="d0186b53-da99-4a0c-a9f3-e686aef25b30" /> creationTime value="20090818163718" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> </pre>

Message Type	Message
	<pre> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="101" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>SIMPSON</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> </patient> </registrationEvent> </subject> </controlActProcess> </subject> </controlActProcess> </pre>

Message Type	Message
	<pre> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>AMSTRONG</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
PIX Query	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="6d79fe84-133c-4c15-9534-79642f13881b" /> creationTime value="20090821104312" /> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> </pre>

Message Type	Message
	<pre> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10502Q1" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> patientIdentifier> value extension="109" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </PRPA_IN201309UV02> </pre>
PIX Query Response Validation Context	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:typeCod e/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- acknowledgementDetail --> </pre>

Message Type	Message
	<pre> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:acknowledgementDetail/@typeCode</XPath> Value> PlainText>E</PlainText> </Value> </DataValueLocationItem> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:acknowledgementDetail/*:code</XPath> Value> PlainText>240</PlainText> </Value> </DataValueLocationItem> <!-- acknowledgementDetail text --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:acknowledgementDetail/*:text</XPath> Value> PlainText>Unknown key identifier</PlainText> </Value> </DataValueLocationItem> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck/*:queryResponseCode/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>

2.2.3 TEST CASE 10503V3: PIX QUERY CASE 4

Message Type	Message
PIX Query	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="5d9728c6-cdb9-4039-a2a5-379203f08099" /> creationTime value="20090821104335" /> </pre>

Message Type	Message
	<pre> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10503Q1" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> dataSource> value root="1.2.3.4.5.8888" /> semanticsText>DataSource.id</semanticsText> </dataSource> patientIdentifier> value extension="100" assigningAuthorityName="UnknownDomain1" root="1.2.3.4.5.9999" /> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </pre>

Message Type	Message
	</PRPA_IN201309UV02>
PIX Query Response Validation Context	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/ *:acknowledgement/ *:typeCod e/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- acknowledgementDetail code --> <!-- acknowledgementDetail --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/ *:acknowledgement/ *:acknowl edgementDetail/@typeCode</XPath> Value> PlainText>E</PlainText> </Value> </DataValueLocationItem> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/ *:acknowledgement/ *:acknowl edgementDetail/ *:code</XPath> Value> PlainText>240</PlainText> </Value> </DataValueLocationItem> <!-- acknowledgementDetail text --> DataValueLocationItem> </pre>

Message Type	Message
	<pre> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:acknowledgementDetail/*:text</XPath> Value> PlainText>Unknown key identifier</PlainText> </Value> </DataValueLocationItem> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck/*:queryResponseCode/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>

2.2.4 TEST CASE 10506V3: PIX QUERY PATIENT UPDATE

Message Type	Message
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/message"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> DataValueLocationItem> <XPath>/*:MCCI_IN000002UV01/*:acknowledgement/*:typeCode/@code</XPath> Value> PlainText>CA</PlainText> </Value> </DataValueLocationItem> </pre>

Message Type	Message
	<pre> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
Patient Add: Number 1	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" > id root="4c117afc-6473-4689-82ec-241e1196a57b" /> creationTime value="20090831133241" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="100" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>Simpson</family> </pre>

Message Type	Message
	<pre> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Albert</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON" /> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Patient Add: Number 2	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="5c2951ed-5896-487b-97e1-16a8031afb46" /> creationTime value="20090831133409" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" </pre>

Message Type	Message
	<pre> classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="200" assigningAuthorityName="domain2" root="1.2.3.4.5.2000" /> statusCode code="active" /> patientPerson> name> given>Jhn</given> family>Spson</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840715" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Albrt</family> </name> </relationshipHolder1> </pre>

Message Type	Message
	<pre> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON" /> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Patient Update	<pre> PRPA_IN201302UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="36834a1b-6d39-456c-8040-6d9a3a25bee6" /> creationTime value="20090831133457" /> interactionId extension="PRPA_IN201302UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </pre>

Message Type	Message
	<pre> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id assigningAuthorityName="domain2" extension="200" root="1.2.3.4.5.2000" /> statusCode code="active" /> patientPerson> name> given>John</given> family>Simpson</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Albert</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON" /> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </pre>

Message Type	Message
	</PRPA_IN201302UV02>
PIX Query	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="28bbddad-5a16-4cf6-a57e-a4579dcc2de6" /> creationTime value="20090831152446" /> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10506Q1" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> dataSource> value extension="200" assigningAuthorityName="domain2" root="1.2.3.4.5.2000" /> semanticsText>DataSource.id</semanticsText> </dataSource> patientIdentifier> value extension="100" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> </pre>

Message Type	Message
	<pre> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </PRPA_IN201309UV02> </pre>
PIX Query Response Validation Context	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- assigningAuthorityName --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@assigningAuthorit yName</XPath> Value> PlainText>domain1</PlainText> </Value> </DataValueLocationItem> <!-- extension --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@extension</XPath> > Value> PlainText>100</PlainText> </Value> </DataValueLocationItem> <!-- root --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:subject/ *:registrationEvent/*:subject1/*:patient/*:id/@root</XPath> Value> </pre>

Message Type	Message
	<pre> PlainText>1.2.3.4.5.1000</PlainText> </Value> </DataValueLocationItem> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:typeCode/@code</XPath> Value> PlainText>AA</PlainText> </Value> </DataValueLocationItem> <!-- Query ack code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck/*:queryId/@root</XPath> Value> PlainText>10506Q1</PlainText> </Value> </DataValueLocationItem> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck/*:queryResponseCode/@code</XPath> Value> PlainText>OK</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>

2.2.5 TEST CASE 10512V3: PIX PATIENT FEED

Message Type	Message
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> </pre>

Message Type	Message
	<pre> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> DataValueLocationItem> <XPath>/*:MCCI_IN000002UV01/*:acknowledgement/*:typeCode /@code</XPath> Value> PlainText>CA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
Add Patient: Number 1	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="29fc4647-c7a3-4520-ba1f-99cd20e17adf" /> creationTime value="20090818164138" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </pre>

Message Type	Message
	<pre> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="120" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>ELLIE</given> family>EPSILON</family> </name> telecom use="HP" value="tel:314.555.1234" /> telecom use="HP" value="tel:314.555.4444" /> administrativeGenderCode code="F" /> birthTime value="19380224" /> addr> streetAddressLine>1 PINETREE</streetAddressLine> city>WEBSTER</city> state>MO</state> postalCode>63119</postalCode> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Barnes</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </pre>

Message Type	Message
	<pre> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Add Patient: Number 2	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="4bcc595c-c336-4fd7-a2ce-70f33721db30" /> creationTime value="20090818164145" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> </pre>

Message Type	Message
	<pre> id extension="120" assigningAuthorityName="domain1" /> statusCode code="active" /> patientPerson> name> given>ELLIE</given> family>EPSILON</family> </name> telecom use="HP" value="tel:314.555.1234" /> telecom use="HP" value="tel:314.555.4444" /> administrativeGenderCode code="F" /> birthTime value="19380224" /> addr> streetAddressLine>1 PINETREE</streetAddressLine> city>WEBSTER</city> state>MO</state> postalCode>63119</postalCode> </addr> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Add Patient: Numbwe 3	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="2499fee7-5469-49d6-8806-f04e8fe9ea18" /> creationTime value="20090818164154" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> </pre>

Message Type	Message
	<pre> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="120" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>ELLIE</given> family>EPSILON</family> </name> telecom use="HP" value="tel:314.555.1234" /> telecom use="HP" value="tel:314.555.4444" /> administrativeGenderCode code="F" /> birthTime value="19380224" /> addr> streetAddressLine>1 PINETREE</streetAddressLine> city>WEBSTER</city> state>MO</state> postalCode>63119</postalCode> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> </pre>

Message Type	Message
	<pre> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Barnes</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>

2.2.6 TEST CASE 10515V3: PATIENT IDENTITY MERGE

Message Type	Message
Acknowledge Validation Content: All	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="VERSION" Level="ERROR" /> MessageFailure Type="LENGTH" Level="IGNORE" /> MessageFailure Type="DATATYPE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE" Level="ERROR" /> MessageFailure Type="MESSAGE_STRUCTURE_ID" Level="ERROR" /> MessageFailure Type="DATA" Level="ERROR" /> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> DataValueLocationItem> <XPath>/*:MCCI_IN000002UV01/*:acknowledgement/*:typeCode </pre>

Message Type	Message
	<pre> /@code</XPath> Value> PlainText>CA</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>
Add Patient: Number 1	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="02273d23-cb98-4629-9055-c4d9bb518267" /> creationTime value="20090909174947" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="101" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> </patient> </subject1> </registrationEvent> </subject> </controlActProcess> </interactionId> </PRPA_IN201301UV02> </pre>

Message Type	Message
	<pre> statusCode code="active" /> patientPerson> name> given>John</given> family>SIMPSON</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Albert</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON" /> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Add Patient: Number 2	<pre> PRPA_IN201301UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"> id root="007dea5c-5933-4098-b180-b6696e7b05eb" /> creationTime value="20090818164650" /> interactionId extension="PRPA_IN201301UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> </pre>

Message Type	Message
	<pre> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SUBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="105" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name> given>Jhn</given> family>Spson</family> </name> telecom use="HP" value="tel:+1-514-555-1234" /> administrativeGenderCode code="F" /> birthTime value="19840711" /> addr> city>Montreal</city> </addr> personalRelationship classCode="PRS"> code codeSystem="2.16.840.1.113883.5.111" codeSystemName="PersonalRelationshipRoleType" code="MTH" displayName="Mother" /> </pre>

Message Type	Message
	<pre> relationshipHolder1 classCode="PSN" determinerCode="INSTANCE"> name xsi:type="urn:PN" xmlns:urn="urn:hl7-org:v3"> family>Albert</family> </name> </relationshipHolder1> </personalRelationship> </patientPerson> providerOrganization classCode="ORG" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> contactParty classCode="CON"> telecom value="SENDER TELECOM GOES HERE!" /> </contactParty> </providerOrganization> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="1.2.840.114350.1.13.99998.8734" /> </assignedEntity> </custodian> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201301UV02> </pre>
Patient Merge	<pre> PRPA_IN201304UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="e05b87c8-3270-4aa1-afd2-56cc0a37cd58" /> creationTime value="20090817115101" /> interactionId root="2.16.840.1.113883.1.6" extension="PRPA_IN201304UV02" /> processingCode code="P" /> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> </pre>

Message Type	Message
	<pre> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> subject typeCode="SBJ"> registrationEvent classCode="REG" moodCode="EVN"> id nullFlavor="NA" /> statusCode code="active" /> subject1 typeCode="SBJ"> patient classCode="PAT"> id extension="101" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> statusCode code="active" /> patientPerson> name /> </patientPerson> </patient> </subject1> custodian typeCode="CST"> assignedEntity classCode="ASSIGNED"> id root="SENDER OID GOES HERE!" /> </assignedEntity> </custodian> replacementOf typeCode="RPLC"> priorRegistration classCode="REG" moodCode="EVN"> subject1 typeCode="SBJ"> priorRegisteredRole classCode="PAT"> id extension="105" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> </priorRegisteredRole> </subject1> </priorRegistration> </replacementOf> </registrationEvent> </subject> </controlActProcess> </PRPA_IN201304UV02> </pre>
Patient Query	<pre> PRPA_IN201309UV02 ITSVersion="XML_1.0" xmlns="urn:hl7-org:v3"> id root="0c6835b6-3587-461a-9273-f0f70d21de4f" /> creationTime value="20090821104419" /> interactionId extension="PRPA_IN201309UV02" root="2.16.840.1.113883.1.6" /> processingCode code="P" /> </pre>

Message Type	Message
	<pre> processingModeCode code="R" /> acceptAckCode code="AL" /> receiver typeCode="RCV"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99999.4567" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99999.1234" /> </representedOrganization> </asAgent> </device> </receiver> sender typeCode="SND"> device classCode="DEV" determinerCode="INSTANCE"> id root="1.2.840.114350.1.13.99998.8734" /> asAgent classCode="AGNT"> representedOrganization determinerCode="INSTANCE" classCode="ORG"> id root="1.2.840.114350.1.13.99998" /> </representedOrganization> </asAgent> </device> </sender> controlActProcess classCode="CACT" moodCode="EVN"> code code="PRPA_IN201309UV02" codeSystem="2.16.840.1.113883.1.6" /> queryByParameter> queryId root="10515Q1" /> statusCode code="new" /> responsePriorityCode code="I" /> parameterList> patientIdentifier> value extension="105" assigningAuthorityName="domain1" root="1.2.3.4.5.1000" /> semanticsText>Patient.id</semanticsText> </patientIdentifier> </parameterList> </queryByParameter> </controlActProcess> </PRPA_IN201309UV02> </pre>
Patient Query Response Validation Context	<pre> <?xml version="1.0" encoding="UTF-8" ?> n1:HL7V3MessageValidationContextDefinition xmlns:n1="http://www.nist.gov/hl7/core/xmlbeans/validation/m essage"> FailureInterpretation> MessageFailure Type="DATA" Level="ERROR" /> </pre>

Message Type	Message
	<pre> MessageFailure Type="MESSAGE_VALIDATION_CONTEXT" Level="USER" /> </FailureInterpretation> MessageInstanceSpecificValues> <!-- ack type code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:acknowledgement/*:typeCod e/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> <!-- Query ack code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryId/@root</XPath> Value> PlainText>10515Q1</PlainText> </Value> </DataValueLocationItem> --> <!-- Query response code --> DataValueLocationItem> <XPath>/*:PRPA_IN201310UV02/*:controlActProcess/*:queryAck /*:queryResponseCode/@code</XPath> Value> PlainText>AE</PlainText> </Value> </DataValueLocationItem> </MessageInstanceSpecificValues> </n1:HL7V3MessageValidationContextDefinition> </pre>